

## Anti-BMP2 antibody (45-60) {Biotin} (STJ500273)

STJ500273

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-BMP2 (45-60) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IF, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Pig

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.729 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:, 1000 DB: 1:2, 000 ELISA: 1:2, 000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">12156</a>
<b>Gene Symbol</b>	<a href="#">Bmp2</a>
<b>Uniprot ID</b>	<a href="#">BMP2_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide corresponding to Mouse BMP2 amino acid region 45-60 conjugated to Keyhole Limpet Haemocyanin (KLH).
<b>Immunogen Region</b>	45-60
<b>Specificity</b>	
<b>Immunogen Sequence</b>	