

Anti-BMPR1A antibody (460-510) {Biotin} (STJ500280)

STJ500280

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BMPR1A (460-510) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	657
Gene Symbol	BMPR1A
Uniprot ID	BMR1A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 460-510 on human Bone morphogenetic protein receptor type-1A.
Immunogen Region	460-510
Specificity	
Immunogen Sequence	