

Anti-BTBD3 antibody (410-460) {FITC} (STJ500300)

STJ500300

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BTBD3 (410-460) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	22903
Gene Symbol	BTBD3
Uniprot ID	BTBD3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 410-460 on human BTBD3 protein.
Immunogen Region	410-460
Specificity	
Immunogen Sequence	