

Anti-Beta crystallin S antibody (120-178) (STJ500307)

STJ500307

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Beta crystallin S (120-178) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	12970
Gene Symbol	Crygs
Uniprot ID	CRYGS_MOUSE
Immunogen	C-epitope Synthetic peptide taken within amino acid region 120-178 on mouse BXTIS protein.
Immunogen Region	120-178
Specificity	
Immunogen Sequence	



Immunohistochemistry of mouse tissue with Anti-Beta crystallin S primary antibody (STJ500307) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.