

Anti-C2orf88 antibody (1-50) {FITC} (STJ500319)

STJ500319

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-C2orf88 (1-50) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:15, 000-1:20, 000 ELISA: 1:15, 000-1:20, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84281
Gene Symbol	C2orf88
Uniprot ID	SMAKA_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 1-50 on human C2orf88 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	