

Anti-Cannabinoid Receptor II antibody (300-350) {Biotin} (STJ500370)

STJ500370

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cannabinoid Receptor II (300-350) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.61-0.63 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10,000 ELISA: 1:10,000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	12802
Gene Symbol	Cnr2
Uniprot ID	CNR2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 300-350 on mouse CB2 Receptor protein common to both CB2L and CB2S.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	