

Anti-CCAR1 antibody (645-695) (STJ500384)

STJ500384

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CCAR1 (645-695) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55749
Gene Symbol	CCAR1
Uniprot ID	CCAR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 645-695 on human Cell division cycle and apoptosis regulator protein 1.
Immunogen Region	645-695
Specificity	
Immunogen Sequence	



Immunohistochemistry of mouse kidney with Anti-CCAR1 primary antibody (STJ500384) at 1:100 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.