

Anti-MAP3K12 antibody (STJ500390)

STJ500390

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MAP3K12 is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence and Western Blot research applications.
Applications	ELISA, ICC, IF, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IF: 1:500
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7786
Gene Symbol	MAP3K12
Uniprot ID	M3K12_HUMAN
Immunogen	Purified Protein tyrosine kinase 7 protein was emulsified in complete adjuvant and used to immunize the rabbits.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



with Anti-MAP3K12 Antibody (STJ500390) of Mouse Embryonic Stem Cells at 1:100 dilution in buffer.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081