

## Anti-MAP3K12 antibody {Biotin} (STJ500391)

STJ500391

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-MAP3K12 is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence and Western Blot research applications.
<b>Applications</b>	ELISA, ICC, IF, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.59 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IF: 1:100
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">7786</a>
<b>Gene Symbol</b>	<a href="#">MAP3K12</a>
<b>Uniprot ID</b>	<a href="#">M3K12_HUMAN</a>
<b>Immunogen</b>	Purified Protein tyrosine kinase 7 protein was emulsified in complete adjuvant and used to immunize the rabbits.
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	