

Anti-MAP3K12 antibody {FITC} (STJ500392)

STJ500392

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MAP3K12 is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence and Western Blot research applications.
Applications	ELISA, ICC, IF, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IF: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7786
Gene Symbol	MAP3K12
Uniprot ID	M3K12_HUMAN
Immunogen	Purified Protein tyrosine kinase 7 protein was emulsified in complete adjuvant and used to immunize the rabbits.
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	