

Anti-CCR4 antibody (1-50) {Biotin} (STJ500406)

STJ500406

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CCR4 (1-50) is suitable for use in Confocal Microscopy, Dot Blot, ELISA, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	CF, DB, ELISA, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Cat, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1233
Gene Symbol	CCR4
Uniprot ID	CCR4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human C-C chemokine receptor type 4 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	