

Anti-CCR6 antibody (1-50) {Biotin} (STJ500416)

STJ500416

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CCR6 (1-50) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.68 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1235
Gene Symbol	CCR6
Uniprot ID	CCR6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human C-C chemokine receptor type 6 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	