

Anti-CD300LF antibody (230-290) {Biotin} (STJ500461)

STJ500461

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CD300LF (230-290) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146722
Gene Symbol	CD300LF
Uniprot ID	CLM1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 230-290 on human CMRF35-like molecule 1/CD300LF.
Immunogen Region	230-290
Specificity	
Immunogen Sequence	