

## Anti-CD300LF antibody (230-290) {FITC} (STJ500462)

STJ500462

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-CD300LF (230-290) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">146722</a>
<b>Gene Symbol</b>	<a href="#">CD300LF</a>
<b>Uniprot ID</b>	<a href="#">CLM1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 230-290 on human CMRF35-like molecule 1/CD300LF.
<b>Immunogen Region</b>	230-290
<b>Specificity</b>	
<b>Immunogen Sequence</b>	