

Anti-CD3G antibody (131-181) {Biotin} (STJ500466)

STJ500466

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CD3G (131-181) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	917
Gene Symbol	CD3G
Uniprot ID	CD3G_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 131-181 on human T-cell surface glycoprotein CD3 gamma protein.
Immunogen Region	131-181
Specificity	
Immunogen Sequence	