

Anti-CD86 antibody (250-329) (STJ500486)

STJ500486

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CD86 (250-329) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	942
Gene Symbol	CD86
Uniprot ID	CD86_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-329 on human CD86 protein.
Immunogen Region	250-329
Specificity	
Immunogen Sequence	



Western blot with Anti-CD86 Antibody (STJ500486) with peptide on at 1:500 dilution in buffer.