

Anti-CGAS antibody (350-450) {Biotin} (STJ500509)

STJ500509

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CGAS (350-450) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.60 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115004
Gene Symbol	CGAS
Uniprot ID	CGAS_HUMAN
Immunogen	Synthetic peptide taken within amino acid regions 350-450 on human cGAS protein.
Immunogen Region	350-450
Specificity	
Immunogen Sequence	