

Anti-CGB1 antibody (142-165) (STJ500518)

STJ500518

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CGB1 (142-165) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.63 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114335
Gene Symbol	CGB1
Uniprot ID	CGB1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 142-165 on human Chorionic gonadotropin beta protein.
Immunogen Region	142-165
Specificity	This antibody will detect CGB1, CGB2, CGB5, CGB7, CGB8.
Immunogen Sequence	