

Anti-CHN2 antibody (1-100) {Biotin} (STJ500532)

STJ500532

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CHN2 (1-100) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68-0.71 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 0 IHC: 0 ICC: 0 CFM: 0
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1124
Gene Symbol	CHN2
Uniprot ID	CHIO_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-100 on human beta-chimaerin isoform 3 protein.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	