

Anti-SHISA9 antibody (250-300) {FITC} (STJ500557)

STJ500557

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SHISA9 (250-300) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64-0.72 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	72555
Gene Symbol	Shisa9
Uniprot ID	SHSA9_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 250-300 on CKAMP44a protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	