

Anti-Calcyon antibody (1-50) (STJ500576)

STJ500576

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Calcyon (1-50) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications. |
| Applications | ELISA, IHC, IP, WB |
| Host/Source | Rabbit |
| Reactivity | Mouse, Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.55-0.75 µg/µl |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 |
| Formulation | |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 192349 |
| Gene Symbol | Caly |
| Uniprot ID | CALY_RAT |
| Immunogen | Synthetic peptide taken within amino acid region 1-50 on rat Calcyon protein. |
| Immunogen Region | 1-50 |
| Specificity | |
| Immunogen Sequence | |



Western blot with Anti-Calcyon Antibody (STJ500576) of mouse brain at 1:500 dilution in buffer.