

Anti-CALCA antibody (30-80) {Biotin} (STJ500587)

STJ500587

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CALCA (30-80) is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Pig, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution	WB: 1:500
Range	ELISA: 1:10, 000 IP: 1:150 IHC: 1:100 ICC: 1:100 IF: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	796
Gene Symbol	CALCA
Uniprot ID	CALC_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 30-80 on human Calcitonin protein.
Immunogen Region	30-80
Specificity	
Immunogen Sequence	