

Anti-SLC31A1 antibody (1-100) {FITC} (STJ500633)

STJ500633

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC31A1 (1-100) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1317
Gene Symbol	SLC31A1
Uniprot ID	COPT1_HUMAN
Immunogen	N-epitope synthetic peptide taken within amino acid region 1-100 on human High affinity copper uptake protein 1.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	