

## Anti-CX3CL1 antibody (350-395) {FITC} (STJ500638)

STJ500638

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-CX3CL1 (350-395) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 µg/µl
<b>Conjugation</b>	FITC
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	20312
<b>Gene Symbol</b>	Cx3cl1
<b>Uniprot ID</b>	X3CL1_MOUSE
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 350-395 on mouse Fractalkine protein.
<b>Immunogen Region</b>	350-395
<b>Specificity</b>	
<b>Immunogen Sequence</b>	