

Anti-CXCR1 antibody (272-350) (STJ500640)

STJ500640

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CXCR1 (272-350) is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, ICC, IF, IHC-P, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	5-.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution	WB: 1:500
Range	DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:250 IHC: 1:100-1:1000 IF: 1:100-1:1000
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3577
Gene Symbol	CXCR1
Uniprot ID	CXCR1_HUMAN
Immunogen	Synthetic peptide, corresponding to a region within C terminal amino acids 272-350 of Human CXCR1.
Immunogen Region	272-350
Specificity	
Immunogen Sequence	