

## Anti-CXCR5 antibody (320-370) {Biotin} (STJ500664)

STJ500664

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-CXCR5 (320-370) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	CM, ELISA, ICC, IF, IHC, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.50-1.20 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	643
<b>Gene Symbol</b>	CXCR5
<b>Uniprot ID</b>	CXCR5_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 320-370 on human Chemokine (C-X-C motif) receptor 5 protein.
<b>Immunogen Region</b>	320-370
<b>Specificity</b>	
<b>Immunogen Sequence</b>	