

Anti-ACKR3 antibody (STJ500669)

STJ500669

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ACKR3 is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.68 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57007
Gene Symbol	ACKR3
Uniprot ID	ACKR3_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human CXCR7 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot with Anti-ACKR3 Antibody (STJ500669) of mouse cerebellum at 1:250 dilution in buffer.