

Anti-DAPK1 antibody (1370-1420) (STJ500726)

STJ500726

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DAPK1 (1370-1420) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1612
Gene Symbol	DAPK1
Uniprot ID	DAPK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1370-1420 on human Death-associated protein kinase 1 protein.
Immunogen Region	1370-1420
Specificity	
Immunogen Sequence	