

Anti-DAPK2 antibody (300-350) {Biotin} (STJ500730)

STJ500730

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DAPK2 (300-350) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13143
Gene Symbol	Dapk2
Uniprot ID	DAPK2_MOUSE
Immunogen	Synthetic peptide take within amino acid region 300-350 on mouse Death-associated protein kinase 2 protein.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	