

Anti-ZIP Kinase antibody (400-448) {Biotin} (STJ500733)

STJ500733

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZIP Kinase (400-448) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13144
Gene Symbol	Dapk3
Uniprot ID	DAPK3_MOUSE
Immunogen	Synthetic peptide taken within amino acid region on 400-448 on mouse Death-associated protein kinase 3 protein.
Immunogen Region	400-448
Specificity	
Immunogen Sequence	