

Anti-SLC6A3 antibody (550-600) {Biotin} (STJ500736)

STJ500736

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC6A3 (550-600) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:150 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6531
Gene Symbol	SLC6A3
Uniprot ID	SC6A3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 550-600 on SLC6A3 protein.
Immunogen Region	550-600
Specificity	
Immunogen Sequence	