

Anti-DCK antibody (210-260) (STJ500741)

STJ500741

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-DCK (210-260) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Chicken
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10,000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13178
Gene Symbol	Dck
Uniprot ID	DCK_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 210-260 on mouse Deoxycytidine kinase.
Immunogen Region	210-260
Specificity	
Immunogen Sequence	