

## Anti-DCK antibody (210-260) {Biotin} (STJ500742)

STJ500742

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Chicken polyclonal antibody anti-DCK (210-260) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IHC, IP, WB
<b>Host/Source</b>	Chicken
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10,000 IP: 1:200 IHC: 1:50-1:100
<b>Formulation</b>	
<b>Isotype</b>	IgY
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">13178</a>
<b>Gene Symbol</b>	<a href="#">Dck</a>
<b>Uniprot ID</b>	<a href="#">DCK_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 210-260 on mouse Deoxycytidine kinase.
<b>Immunogen Region</b>	210-260
<b>Specificity</b>	
<b>Immunogen Sequence</b>	