

Anti-Deformed Epidermal Autoregulatory Factor 1 antibody (500-565) (STJ500747)

STJ500747

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Deformed Epidermal Autoregulatory Factor 1 (500-565) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:150 ICC: 1:50-1:150
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	54006
Gene Symbol	Deaf1
Uniprot ID	DEAF1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 500-565 on mouse Deformed epidermal autoregulatory factor 1 homolog protein.
Immunogen Region	500-565
Specificity	
Immunogen Sequence	



Immunohistochemistry of rat tissue with Anti-Deformed Epidermal Autoregulatory Factor 1 primary antibody (STJ500747) at 1:50 dilution in Buffer. Section was treated with DAB and Haematoxylin stain. Magnification is at 40X.