

Anti-DFos antibody (5-80) (STJ500756)

STJ500756

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DFos (5-80) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IP, WB
Host/Source	Rabbit
Reactivity	D.melanogaster

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:1, 000 IP: 1:250 IHC: 1:50-1:200 ICC: 1:50-1:200 IF: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 5-80 on DFos protein.
Immunogen Region	5-80
Specificity	
Immunogen Sequence	



Immunofluorescence of third instar ventral ganglion with Anti-DFos Antibody (STJ500756) at 1:250 dilution in Buffer.