

Anti-DNER antibody (700-737) {Biotin} (STJ500778)

STJ500778

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DNER (700-737) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	92737
Gene Symbol	DNER
Uniprot ID	DNER_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-737 on human Delta and Notch-like Epidermal growth factor-related receptorprotein.
Immunogen Region	700-737
Specificity	
Immunogen Sequence	