

Anti-DRAXIN antibody (300-330) {Biotin} (STJ500790)

STJ500790

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DRAXIN (300-330) is suitable for use in ELISA, Immunocytochemistry and Western Blot research applications.
Applications	ELISA, ICC, WB
Host/Source	Rabbit
Reactivity	Bovine, Chicken, Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.64 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:1, 000 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	374946
Gene Symbol	DRAXIN
Uniprot ID	DRAXL_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-330 on C terminal region of Human Draxin protein. Region: NNKCFDDCMC VEGLRCYAKF HRNRRVTRRK
Immunogen Region	300-330
Specificity	
Immunogen Sequence	