

Anti-DRP1 antibody (450-490) {Biotin} (STJ500793)

STJ500793

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DRP1 (450-490) is suitable for use in ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, ICC, IF, IHC, IP, WB
Host/Source	Rabbit
Reactivity	D.melanogaster, Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.68 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:450 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74006
Gene Symbol	Dnm1l
Uniprot ID	DNNM1L_MOUSE
Immunogen	Synthetic peptide within 450-490 on human DRP1 protein
Immunogen Region	450-490
Specificity	
Immunogen Sequence	