

Anti-DPYSL4 antibody (509-527) (STJ500795)

STJ500795

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DPYSL4 (509-527) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Cat, Human, Monkey, Mouse, Pig, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:450 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10570
Gene Symbol	DPYSL4
Uniprot ID	DPYL4_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 509-527 of Mouse CRMP3. The peptide was covalently modified to achieve the desired antigenicity before coupling to carrier protein.
Immunogen Region	509-527
Specificity	
Immunogen Sequence	



Western blot with Anti-DPYSL4 Antibody (STJ500795) with A) CRMP3 Positive Control B) SKNB21 at 1:500 dilution in buffer.