

Anti-DWORF antibody (Full Length) {Biotin} (STJ500799)

STJ500799

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DWORF (Full Length) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Strit1
Uniprot ID	DWORF_MOUSE
Immunogen	Synthetic peptide corresponding to full length Sarcoplasmic/endoplasmic reticulum calcium ATPase regulator DWORF protein.
Immunogen Region	Full Length
Specificity	
Immunogen Sequence	