

Anti-E2F2 antibody (70-120) {FITC} (STJ500821)

STJ500821

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E2F2 (70-120) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.75-1.25 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	242705
Gene Symbol	E2f2
Uniprot ID	E2F2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 70-120 on E2F2 protein.
Immunogen Region	70-120
Specificity	
Immunogen Sequence	