

Anti-E2F5 antibody (5-55) {Biotin} (STJ500829)

STJ500829

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E2F5 (5-55) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.88-0.96 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13559
Gene Symbol	E2f5
Uniprot ID	E2F5_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 5-55 on E2F5 protein.
Immunogen Region	5-55
Specificity	
Immunogen Sequence	