

Anti-EHD1 antibody (515-565) {Biotin} (STJ500850)

STJ500850

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EHD1 (515-565) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10938
Gene Symbol	EHD1
Uniprot ID	EHD1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 515-565 on human EH-domain containing protein 1.
Immunogen Region	515-565
Specificity	
Immunogen Sequence	