

Anti-EHD2 antibody (420-450) {FITC} (STJ500854)

STJ500854

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EHD2 (420-450) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	30846
Gene Symbol	EHD2
Uniprot ID	EHD2_HUMAN
Immunogen	Synthetic peptide taken within amino acid 420-450 on human EH-Domain Containing Protein 2.
Immunogen Region	420-450
Specificity	
Immunogen Sequence	