

## Anti-CD39 antibody (380-430) (STJ500891)

STJ500891

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-CD39 (380-430) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">12495</a>
<b>Gene Symbol</b>	<a href="#">Entpd1</a>
<b>Uniprot ID</b>	<a href="#">ENTP1_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 380-430 on mouse EctoNucleoside TriPhosphate Diphosphohydrolase 1 protein.
<b>Immunogen Region</b>	380-430
<b>Specificity</b>	
<b>Immunogen Sequence</b>	