

Anti-RAPGEF3 antibody (526-541) (STJ500901)

STJ500901

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RAPGEF3 (526-541) is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IHC, WB
Host/Source	Rabbit
Reactivity	Horse, Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.59 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10411
Gene Symbol	RAPGEF3
Uniprot ID	RPGF3_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 526-541 of Epac1 protein, peptide was covalently modified. Sequence DHSVLTQLQLPVTASVR
Immunogen Region	526-541
Specificity	Does not cross react with Epac2.
Immunogen Sequence	