

Anti-Ephrin Receptor B3 antibody (100-200) {Biotin} (STJ500937)

STJ500937

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ephrin Receptor B3 (100-200) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.58-0.61 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:750 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13845
Gene Symbol	Ephb3
Uniprot ID	EPHB3_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 100-200 on EphBR3 protein.
Immunogen Region	100-200
Specificity	
Immunogen Sequence	