

Anti-EPHB4 antibody (900-984) (STJ500939)

STJ500939

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EPHB4 (900-984) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.62 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2050
Gene Symbol	EPHB4
Uniprot ID	EPHB4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-984 on human Ephrin B4 receptor protein.
Immunogen Region	900-984
Specificity	
Immunogen Sequence	



Western blot with Anti-EPHB4 Antibody (STJ500939) of SKNMC cell at 1:500 dilution in buffer.