

Anti-EPHB4 antibody (900-984) {FITC} (STJ500941)

STJ500941

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EPHB4 (900-984) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.62 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2050
Gene Symbol	EPHB4
Uniprot ID	EPHB4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-984 on human Ephrin B4 receptor protein.
Immunogen Region	900-984
Specificity	
Immunogen Sequence	