

## Anti-ERG antibody (400-450) {Biotin} (STJ500950)

STJ500950

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-ERG (400-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.55-0.75 µg/µl
<b>Conjugation</b>	Biotin
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 DB: 1:10,000 ELISA: 1:10,000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">2078</a>
<b>Gene Symbol</b>	<a href="#">ERG</a>
<b>Uniprot ID</b>	<a href="#">ERG_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 400-450 on human Transcriptional regulator ERG protein.
<b>Immunogen Region</b>	400-450
<b>Specificity</b>	This is a Pan ERG antibody and will label isoforms 1, 2, 3, and 4.
<b>Immunogen Sequence</b>	